# Blade<sup>XT</sup>

## Instructions for Use

BXTAG1SF-BXTAG9SF  
BXTAG1SM-BXTAG9SM

EN	Instructions for Use	2
DA	Brugsanvisning	13
NO	Bruksanvisning	24
FI	Käyttöohjeet	35
SV	Bruksanvisning	46
RU	Инструкция протезиста	57
ZH	使用说明	68
AR	إرشادات الاستخدام	79

# Contents

Contents .....	2
1 Description and Intended Purpose .....	3
2 Safety Information .....	4
3 Construction .....	5
4 Alignment .....	6
4.1 Static Alignment .....	6
4.2 Dynamic Alignment .....	6
5 Fitting Advice .....	7
6 Maintenance .....	7
7 Limitations on Use .....	8
8 Heel Spring Removal & Replacement .....	9
9 Technical Data .....	10
10 Ordering Information .....	11

---

# 1 Description and Intended Purpose

These Instructions for Use are intended for use by the practitioner and user unless otherwise stated.

The term *device* is used throughout these instructions for use to refer to BladeXT.

Please read and ensure you understand all instructions for use, in particular all safety information and maintenance instructions.

## Application

This device is to be used exclusively as part of a lower limb prosthesis.

Intended for a single user.

This device is designed specifically for general sport and running.

## Activity Level

This device is recommended for users that have the potential to achieve Activity Level 4 (weight restrictions apply, see Section 9 *Technical Data*). Of course there are exceptions and in our recommendation we want to allow for unique, individual circumstances and any such decision should be made with sound and thorough justification.

### Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited household ambulator.

### Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low-level environmental barriers such as curbs, stairs, or uneven surfaces. Typical of the limited community ambulator.

### Activity Level 3

Has the ability or potential for ambulation with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

### Activity Level 4

Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

## Contraindications

This device may not be suitable for Activity Level 1 and 2 users.

## Clinical Benefits

- Allows variable running speeds
- Sole plate improves perceived stability, compared to other running-specific prostheses
- Users demonstrate confidence in prosthetic loading during high activity

## Spring Set Selection

User	User Weight									kg lb	
	44-52 100-115	53-59 116-130	60-68 131-150	69-77 151-170	78-88 171-195	89-100 196-220	101-116 221-255	117-130 256-285	131-147 286-325		148-166 326-365
Jogger	1	2	3	4	5	6	7	8	9	9	Toe and Heel spring set
Runner	2	3	4	5	6	7	8	9	9		
Jogger	Occasional to frequent recreational running										
Runner	Frequent distance running and endurance training										

**Important:** For higher impact users, do not exceed the weight limit for individual springs.

## 2 Safety Information

This warning symbol highlights important safety information which must be followed carefully.

Any changes in the performance or function of the limb e.g. restricted movement, non-smooth motion or unusual noises should be immediately reported to your service provider.

Always use a hand rail when descending stairs and at any other time if available.

To minimize the risk of 'over-stepping' when descending stairs/steps, be aware of the 'C'-shaped toe spring extending behind you. Take care to avoid clipping the edge or back of a stair/step with the spring during descent which could cause you to trip or fall.

Assembly, maintenance and repair of the device must only be carried out by a suitably qualified practitioner.

The device is designed for prolonged submersion and suitable for immersion in fresh water only. Ensure any use of the device in water complies with the conditions given in Section 7 *Limitations on Use*.

The user should contact their practitioner if their condition changes.

Ensure only suitably retrofitted vehicles are used when driving. All persons are required to observe their respective driving laws when operating motor vehicles.

To minimize the risk of slipping and tripping, the sole cover should not be excessively worn, torn or loose.

Be aware of finger trap hazard at all times.


Check sole spring bolts for tightness before use.

# 3 Construction

## Principal Parts:


- Heel & Toe Springs (e-Carbon)
- Heel attachment Screws & Washers (St. Stl. - Black)
- Sole Cover & Wedge (PU)

### Adapter/Male BXTAG1SM-BXTAG9SM

938344


35 Nm  
5/16 Whitworth

### Adapter/ Female BXTAG1SF-BXTAG9SF

938343

35 Nm  
5/16 Whitworth

## 4 Alignment

*The instructions in this section are for practitioner use only.*


### 4.1 Static Alignment

#### Transfemoral Alignment

Align trans-femoral components according to fitting instructions supplied with the knee. Keep the build line relative to this device as shown.

#### Setup Length

This device should be set up with a 10mm heel wedge so that its top surface is horizontal. The length of this device should be 10–15 mm longer than the user's everyday prosthesis with running trainer footwear to allow for spring compression.

#### Build Line

The build line should fall in the middle of the slot on top of this device. This is defined by a vertical alignment indicator line on each side.

### 4.2 Dynamic Alignment

#### Coronal Plane

Ensure that M-L thrust is minimal by adjusting relative positions of socket and foot.

#### Sagittal Plane

On standing the toe should deflect so that the heel is in contact with the ground. On walking check for a smooth transition from heel strike to toe-off.

#### 1. Optimize plantar flexion and dorsiflexion

On running, a heel strike should always be present but this can be very small and subtle. To optimize running gait, plantarflex the foot incrementally in small amounts at a time ensuring a heel strike is maintained. If the foot is too plantarflexed the sole will 'slap' against the ground. To prevent 'slap' dorsiflex the foot.

#### 2. Stiffness and A-P shift

When an optimal plantar flexion/ dorsiflexion angle has been achieved, if when running the toe feels too soft, then shift the weight line posteriorly on the foot. If the toe feels too stiff, then shift the weight line anteriorly on the foot.

#### 3. If the user is sinking at heel strike the heel is too soft

Fitting the wedge supplied will have the effect of stiffening the heel spring, this can be taped in place for trial (see diagram). For permanent fitting the wedge should be glued in place by application of Loctite 424 (926104) on the upper surface of the wedge.

Ensure the top surface of the blade is horizontal to allow for compression in the main Toe spring when loaded with the user's body weight.

---

## 5 Fitting Advice

*The instructions in this section are for practitioner use only.*

Springs are supplied as matched sets i.e. the shin/toe and heel springs are designed to work together to give smooth progression for most users.

Condition	Symptoms	Remedy
Heel too soft	Sinking at heel strike	<ol style="list-style-type: none"><li>1. Add heel wedge</li><li>2. Slide socket anteriorly</li></ol>
Heel too hard	Rapid transition from heel strike through stance phase Difficulty in controlling heel action, foot jars into mid-stance	<ol style="list-style-type: none"><li>1. Remove heel wedge (if fitted)</li><li>2. Slide socket posteriorly</li><li>3. Order a softer Heel Spring category</li></ol>
Toe too soft	Rapid progression through mid-stance Drop-off during running	<ol style="list-style-type: none"><li>1. Move socket posteriorly in relation to foot</li><li>2. Plantarflex foot slightly - note: heel strike must be maintained</li></ol>
Toe too hard	Feels like climbing over toe	<ol style="list-style-type: none"><li>1. Move socket anteriorly in relation to foot</li><li>2. Dorsiflex foot slightly</li></ol>

Please contact your supplier if it is not possible to achieve a smooth gait after following the recommendations above.

---

## 6 Maintenance

Visually inspect the device before any activity, especially the condition of the sole and toe spring, and the tightness of the sole spring bolts.

Excessive wear of the sole cover could lead to loss of grip. Do not use if the sole cover is excessively worn, torn or loose.

Excessive wear or damage to the spring elements could lead to structural failure.

Report any changes in performance of this device to the practitioner/service provider e.g. reduced energy return, unusual noises, softening of the toe and/or heel, or significant wear.

Inform the practitioner/service provider of any changes in body weight and/or activity level.

### Cleaning

Use a damp cloth and mild soap to clean the outside surfaces. DO NOT use aggressive cleansers. Thoroughly dry before use.

*The remaining instructions in this section are for practitioner use only.*

Advise the user that the device should **not** be used if the sole cover is excessively worn, torn or loose. The user should also be advised to visually inspect this device especially the condition of the heel and toe springs before any activity.

---

## 6 Maintenance (cont.)

*The instructions in this section are for practitioner use only.*

It is recommended that the following maintenance is carried out annually:

- Check the sole cover for damage or wear and replace the Heel Spring assembly if necessary.
- Check all bolts for tightness, clean and reassemble as necessary.
- Visually check both the Heel and Toe blade springs for signs of delamination or wear and replace if necessary. Some minor surface damage may occur after a period of use, this does not affect the function or strength of the foot.

The user should be advised of the following:

- Beware of slip hazard at all times especially on wet/smooth surfaces.
- Any changes in performance of this device must be reported to the practitioner e.g. reduced energy return, unusual noises or softening of the toe or heel.
- The practitioner must also be informed of any changes in body weight and/or activity level.
- Excessive wear of the sole cover could lead to loss of grip.
- Excessive wear or damage to the spring elements could lead to structural failure.

---

## 7 Limitations on Use

### Intended Life

A local risk assessment should be carried out based upon activity and usage.

### Lifting Loads

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

### Environment

This device is waterproof to a maximum depth of 1 meter.

Thoroughly rinse this device with fresh water after use in abrasive environments such as those that contain sand or grit, for example, to prevent wear or damage.

Thoroughly rinse with fresh water after use in salt or chlorinated water.

Exclusively for use between -15 °C and 50 °C (5 °F to 122 °F).


Suitable for submersion

## 8 Heel Spring Removal & Replacement


*The instructions in this section are for practitioner use only.*

## 9 Technical Data

Material:	Aluminum, St. Stl. Titanium, Glass Fiber/Carbon Fiber
Operating and Storage Temperature Range:	-15 °C to 50 °C (5 °F to 122 °F)
Component Weight:	850 g (1 lb 14 oz)
Activity Level:	4
Maximum User Weight:	166 kg (366 lb)
Proximal Alignment Attachment:	Male pyramid adapter (Blatchford) Female pyramid adapter (Blatchford)
Build Height:	Male pyramid adapter 230 mm Female pyramid adapter 250 mm
Proximal pyramid to ground level [See diagram below]	

### Fitting length

---

# 10 Ordering Information

## Order Example

<b>BXTAG</b>	<b>3</b>	<b>S</b>	<b>F</b>
	Spring Set (See table below.)		Adapter Type (M/F)

Specify Spring Set from 1 to 9 and adapter type, female [F] or male [M].

BXTAG1SF to BXTAG9SF

BXTAG1SM to BXTAG9SM

e.g. BXTAG3SF, BXTAG6SM

---

## BXT Heel Spring Kit

Includes the overmoulded cover, fasteners and heel wedge.

Spring set	Heel type	Part no.
1, 2	Soft Sports	539083SS
3, 4, 5	Medium Sports	539083MS
6, 7	Firm Sports	539083FS
8, 9	EX Firm Sports	539083EFS

## Proximal adapter

Type	Part no.
Male	189427
Female	189527

## Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

## CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: [www.blatchford.co.uk](http://www.blatchford.co.uk)

Medical Device

Single Patient – multiple use

## Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a practitioner.

## Warranty

This device is warranted for 24 months (excludes paintwork and sole cover).

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions.

See the Blatchford website for the current full warranty statement. Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

## Environmental Aspects

To help prevent potential harm to the environment or to human health from uncontrolled waste disposal, Blatchford offers a take-back service. Please contact Customer Services for details.

## Retaining the Packaging Label

The practitioner is advised to keep the packaging label as a record of the device supplied.

## Trademark Acknowledgements

BladeXT and Blatchford are registered trademarks of Blatchford Products Limited.

## Manufacturer's Registered Address

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**Blatchford Europe GmbH**
**Endolite India Ltd.**