Assembly of Cylinder to Carrier Assembly

Assembly Instructions

932285, 932286,
MH30, MHPYR

| | | |
|----|-------------------------|----|
| EN | Assembly Instructions | 2 |
| DA | Monteringsanvisninger | 6 |
| NO | Monteringsinstruksjoner | 10 |
| FI | Kokoonpano-ohjeet | 14 |
| SV | Monteringsanvisningar | 18 |
| RU | Сборочные инструкции | 22 |
| ZH | 装配说明 | 26 |
| AR | إرشادات التجميع | 30 |

Contents

| | |
|------------------------------------|---|
| Contents | 2 |
| 1 Introduction | 3 |
| 2 Function..... | 3 |
| 3 Lower (Distal) Attachment | 4 |
| 4 Upper (Proximal) Attachment..... | 5 |

1 Introduction

These Assembly Instructions are intended for the practitioner and provide a supplement to the Instructions For Use.


This document details how to correctly fit a Blatchford Hydraulic Control Cylinder to the Mercury range of Carrier Assemblies.

2 Function

The Hydraulic Control Cylinder forms part of the knee frame assembly and fits between the chassis posterior pivot and the carrier distal pivot.

The combined cylinder and carrier assembly is a Swing-and-Stance device, providing:


- Adjustable hydraulic swing phase control to accommodate variable cadence.
- Adjustable yielding hydraulic stance support, activated from knee extension and released by a hyper-extension moment at any time but usually occurring after mid-stance.

3 Lower (Distal) Attachment

First, attach the lower end of the Hydraulic Control Cylinder to the distal end of the Carrier Assembly using the following procedure:

1. Apply Loctite 243 to the Stub Axle threaded holes on the cylinder.
2. Place the Hydraulic Control Cylinder inside the Carrier Assembly, ensuring the piston end is upper-most.
3. Turn the cylinder to align the threaded holes with the holes in the Carrier Assembly.
4. Lubricate the Stub Axles with Lithium Grease (928003) avoiding the threaded end.
5. Insert Stub Axles ensuring selected washers are in place. Select washer size to minimize side-play and ALWAYS use the same color/size washer on both sides. Tighten to the correct torque setting.


| | | | | | |
|--------------------|--|---|--|---|-----------|
| WHITE 0.6 THICK | | + | | ✓ | Preferred |
| BLACK 0.5 THICK | | + | | ✓ | |
| | | + | | ✗ | |

Note... Choose washer thickness to minimise side-play.
After fitting, limb must return to full extension after knee flexion.

4 Upper (Proximal) Attachment

Attach the upper end of the Hydraulic Control Cylinder to the proximal end of the Carrier Assembly using the following procedure:

1. Ensure Trunnion pivot holes in the chassis assembly are not obscured.
2. Apply STP Oil/Mineral Grease (928105) to both washers.
3. Place the washers on either side of the Trunnion and insert between the chassis assembly pivot holes. Align holes, Trunnion and washers.
4. Lubricate the Upper Attachment Pin with Lithium Grease (928003) before pushing into the pivot holes using a Hex Key and a slow twisting action.
5. Flex and extend the knee three or four times to centralize the cylinder.
6. Apply Loctite 222 to the Trunnion Grub Screw and gently tighten. **DO NOT OVERTIGHTEN.**

| | |
|--|---|
| Indhold..... | 6 |
| 1 Introduktion | 7 |
| 2 Funktion..... | 7 |
| 3 Nedre (distal) kobling | 8 |
| 4 Øverste (proksimal) tilkobling | 9 |

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Blatchford Europe GmbH
Endolite India Ltd.