Hydraulic Knee Control

Instructions for Use

932281, 932283, 932285, 932286

| | | |
|----|----------------------------|-----|
| EN | Instructions for Use | 2 |
| FR | Instructions d'utilisation | 18 |
| DE | Gebrauchsanweisung | 34 |
| IT | Istruzioni per l'uso | 50 |
| ES | Instrucciones de uso | 66 |
| NL | Gebruiksaanwijzing | 82 |
| PL | Instrukcja użytkowania | 98 |
| PT | Instruções de utilização | 114 |
| CS | Návod k použití | 130 |

Contents

- Contents2
- 1 Description and Intended Purpose3
- 2 Safety Information4
- 3 Construction5
- 4 Function.....6
- 5 Maintenance6
- 6 Limitations on Use7
- 7 Alignment Procedure.....8
 - 7.1 Static Alignment.....9
 - 7.2 Dynamic Alignment: Swing Resistance..... 11
 - 7.3 Stance Selection Lever/Stirrup Function and Features 13
 - 7.4 Static Alignment Check..... 14
 - 7.5 Dynamic Adjustment..... 15
- 8 Technical Data 16
- 9 Ordering Information 16

1 Description and Intended Purpose

These Instructions for Use are intended for the practitioner and user unless otherwise stated.

The term *device* is used throughout these instructions for use to refer to the Blatchford Hydraulic Knee Control.

Please read and ensure you understand all instructions for use, in particular all safety information and maintenance instructions.

Application

These devices are intended for use with Blatchford Mercury Knee and KX06 (not KX06V2).

Features

- Hydraulic Control Cylinder
- Variable Cadence
- Compact Design
- Adjustable Stance Yield Resistance
- Swing-Only Mode (Suitable for Cycling)
- Cylinder Flexion Lock

Activity Level

This device is suitable for Activity Levels 3 and 4 (weight limits apply, see *Technical Data*). The device is not suitable for Activity Levels 1 and 2 or for use in competitive sports events. These types of users might be better served by a specially designed prosthesis that is optimized for their needs. Of course there are exceptions and in our recommendation we want to allow for unique, individual circumstances and any such decision should be made with sound and thorough justification.

Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited household ambulator.

Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low-level environmental barriers such as curbs, stairs, or uneven surfaces. Typical of the limited community ambulator.

Activity Level 3

Has the ability or potential for ambulation with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

Activity Level 4

Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

Clinical Benefits

- Independently adjustable to suit amputee walking style
- Knee resistance yield allows supported stair descent.
- Knee resistance yield allows supported sitting.
- Flexion lock allows extended standing support.
- Selectable Swing-Only mode to suit some sports activities.

Contraindications

This device is only intended for use as a component of Blatchford Mercury High Activity (MHA) and KX06 (not KX06v2) prosthetic limb products.

2 Safety Information

This warning symbol highlights important safety information which must be followed carefully.

Any changes in the performance or function of the limb e.g. restricted or excessive movement, non-smooth motion or unusual noises should be immediately reported to your service provider.

Always use a hand rail when descending stairs and at any other time if available.

After continuous use the device casing may become hot to the touch.

Any excessive changes in heel height after finalization of alignment may adversely affect limb function.

Do not place near any heat source. Do not leave in direct sunshine or inside a car in hot weather.

The device can be used as shower prosthesis but is not designed for prolonged submersion in water. After contact with water wipe dry immediately. If this device comes into contact with salt or chlorinated water, it should be rinsed with fresh water and dried. Ensure any use of the device complies with the conditions given in *Limitations on Use*.

The device is not suitable for extreme sports, running or cycle racing, ice and snow sports, extreme slopes and steps. Any such activities undertaken are done so completely at the users' own risk. Recreational cycling is acceptable.


Assembly, maintenance and repair of the device must only be carried out by a suitably qualified practitioner.

3 Construction

Principal Parts

- Cylinder Body Aluminium alloy
- Flexion Cap Aluminum alloy
- Piston Rod Stainless steel
- Trunnion Stainless steel
- Dust Cap Polyacetal homopolymer
- Stance Selector Lever/Stirrup Stainless steel
- Bushing Brass (KX06) or aluminum (MHA)
- Internals Aluminum alloy, stainless steel, acetal homopolymer, brass, NBR, hydraulic fluid

Component Identification

4 Function

The device is intended to sit within a frame assembly, fitted between the chassis posterior pivot and the carrier distal pivot.

The combined cylinder and carrier assembly is a Swing-and-Stance device, providing:

- Adjustable hydraulic swing phase control to accommodate variable cadence.
- Adjustable yielding hydraulic stance support, activated from knee extension and released by a hyper-extension moment at any time but usually occurring after mid-stance.

5 Maintenance

No specific maintenance advice is required for this device.

Visually check the device regularly.

The user must be advised:

- Any change in performance or function of the device must be reported to the practitioner e.g. unusual noises*, faster or slower swing or reduced stance support.
- After running or periods of high activity the cylinder may become very warm, this is normal.

Note... * Due to the nature of the hydraulics there may be some slight air noise from the cylinder during the first few steps. This is not detrimental to the function of the unit and should dissipate quickly. If symptoms persist, please consult your practitioner.

Cleaning

Use a damp cloth and mild soap to clean outside surfaces, do not use aggressive cleansers.

Ensure the user has read and understood all safety and user-level maintenance information.

Advise the user that a regular visual check of the device is recommended and signs of wear that may affect function should be reported to their service provider.

Storage & Handling

When storing for prolonged periods place cylinder vertically with trunnion uppermost.

Alternatively and in addition, air management may be assisted by cycling the cylinder several times with the stance resistance switched off.

Use product packaging supplied.

6 Limitations on Use

Intended Life

A local risk assessment should be carried out based upon activity and usage.

Lifting Loads

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

Environment

This device can be used as shower prosthesis but is not designed for prolonged submersion in water. After contact with water wipe dry immediately. If this device comes into contact with salt or chlorinated water, it should be rinsed with fresh water and dried. Light surface corrosion will not affect the function or security of this device. However, if heavy corrosion is evident, stop using this device and contact your practitioner.

Exclusively for use between -10 °C and 50 °C (14 °F to 122 °F).

Please refer to the original KX06 and/or Mercury Knee product IFUs for guidance on outdoor use.

Suitable for shower use

7 Alignment Procedure

The instructions in this section are for practitioner use only.

When aligning the prosthesis for use with this device, positioning of the build line is critical.


Knees should be aligned so the build line passes up to 5 mm behind the knee centre (posterior).

Knees only release from stance after a hyper-extension moment at the knee (which occurs naturally in late stance), therefore flexion should be initiated whilst weight is still on the limb.

The notes above also apply to the original KX06 (not KX06V2) knee (not shown).

Mercury knees may be fitted with a distal pyramid instead of the tube clamp to allow the addition of pyramid-based components. Pyramids should, however, be aligned so the shin components are vertical to maximise the function of the knee.

These notes also apply to the Mercury Hi-Activity Discontinuous shin assembly (not shown).


7.1 Static Alignment

To achieve the optimal function from the device within a Mercury Knee product, the knee must be aligned geometrically unstable. For other configurations please refer to associated documentation.

Check flexion is fully accommodated when worn by patient.

Checking Anterior-Posterior (A-P) Alignment


Adjusting Yield/Stance Resistance

2


Having first checked A-P alignment, see Box 1 (initial adjustment position is factory set):

- a) Ensure lever (stirrup) is in DOWN position
- b) Standing with feet side-by-side, try to sit (do not hyperextend limb)


If excessive resistance is felt, turn the Yield Adjustment Screw anti-clockwise, using the Hex Key, until the correct resistance is felt to sitting naturally.


Stance Resistance (Yield). Do not turn the yield adjustment screw counterclockwise beyond 180° position as this will cause the control unit to operate inconsistently (may or may not release from stance).

Yield Adjustment Screw

9 = Maximum Stance Resistance


1 = Minimum Stance Resistance

7.2 Dynamic Alignment: Swing Resistance

Adjusting Flexion Resistance

3

- Initial adjustments are factory set (Extension 2, Flexion 4)
- Observe the amputee walking
- If there is excessive heel rise, increase resistance
- If there is insufficient heel rise, decrease resistance


The Flexion Cap is numbered from left (MIN/1) to right (MAX/10).

Directly below the Flexion Cap is a fixed Indicator Mark.

When the cap is turned counter-clockwise, when MIN/1 is over the mark, flexion resistance is at minimum.

When the cap is turned clockwise, when MAX/10 is over the mark, flexion resistance is at maximum.

Warning! If adjusting the flexion resistance has no effect on heel rise, check stirrup is in down position and ensure that the amputee is initiating flexion whilst still loading the toe.


Remember a hyperextension moment is needed about the knee to initiate flexion.

Adjusting Extension Resistance

4

- a) Observe the amputee walking
- b) If there is excessive terminal impact on knee extension, increase resistance
- c) If the knee does not extend satisfactorily, decrease resistance

Note... As a 'rule of thumb', flexion resistance should be higher than extension resistance.


To adjust extension resistance use the same Flexion Cap (make a note of the flexion resistance before you start). By turning the cap in the clockwise direction until it reaches maximum and then continuing, the indicator mark will begin moving from right to left.

The indicator mark will move clockwise from 1 to 10 on the extension marker which is located directly below the indicator mark.

The flexion adjuster cap should then be readjusted to the desired flexion resistance.

Note... To reduce extension resistance turn the adjuster cap to minimum and continue turning until the desired extension resistance is reached. Readjust the flexion adjustment cap accordingly.

Minimum
Extension
Resistance

Maximum
Extension
Resistance


Note... If there is any doubt over the position of the indicator ring in relation to the settings, its position can be reset by turning the adjustment cap to maximum flexion, then maximum extension prior to making any adjustments.

Note... If using adjuster tool 940091, take care not to force the adjuster ring and overcome the adjustment stops as this could cause loss of function.

7.3 Stance Selection Lever/Stirrup Function and Features

Note... Always confirm function when changing modes.


How to Switch Off Stance Resistance (for Cycling, etc.)


- Load toe of prosthesis
- Lift lever (stirrup)

 **This will override stance control.**

- Stance now off and knee will be unstable
- Lower lever (stirrup) to resume normal operation


How to Lock Knee Against Flexion (for Prolonged Standing, etc.)


- Flex knee slightly, without hyperextending the knee first
- Lift lever (stirrup)

 **This will override stance control.**


- The knee is now locked against flexion but will extend
- Lower lever (stirrup) to resume normal operation

7.4 Static Alignment Check

7.5 Dynamic Adjustment

8 Technical Data

Operating and Storage
Temperature Range:

-10 °C to 50 °C
(14 °F to 122 °F)

Component Weight:

475 g (1.05 lb)


Activity Level:

3-4

Maximum User Weight:

Limitations apply -
please see carrier details

Key Dimensions (See Diagram):

9 Ordering Information

New and Service Exchange Units

| New | S/Ex | Description | Trunnion |
|--------|--------|--|----------|
| 932281 | 932283 | Swing-and-Stance Cylinder for KX06 | 8 mm |
| 932285 | 932286 | Swing-and-Stance Cylinder for Mercury Knee | 10 mm |

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: www.blatchford.co.uk

Medical Device

Single Patient – multiple use

Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a practitioner.

Warranty

New devices are warranted for 36 months.

Service Exchange devices are warranted for 12 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions.

See the Blatchford website for the current full warranty statement.

Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

Environmental Aspects

This product contains hydraulic oil, mixed metals and plastics. It should be recycled where possible in accordance with local waste recycling regulations.

Retaining the Packaging Label

The Practitioner is advised to keep the packaging label as a record of the device supplied.

Trademark Acknowledgements

Blatchford is a registered trademark of Blatchford Products Limited.

Manufacturer's Registered Address

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Blatchford Europe GmbH
Endolite India Ltd.