

Paediatric Footwear

Issue 4

KSA

 **Mendivil**

Mendivil Paediatric Footwear

Mendivil Paediatric Footwear is exclusively distributed by Blatchford.

It features a range of styles in sizes 16 to 45 with two availability options: Stocked Footwear or Made to Order Footwear. Stocked Footwear is available for a 24 hour delivery and limited to the three styles, colour and fastenings as shown on page 3, which are ideal for school.

Made to Order Footwear styles are shown from page 11 onwards. This range offers a wide selection of colour options in either lace or reverse hook and loop fastening in two back heights: UK and Continental. Standard, anti-varus and orthoses accommodating width fittings plus three stiffener options are also available. Not all options are applicable to every style, please check the specification provided alongside each boot in the catalogue on the relevant page.

All Mendivil Paediatric Footwear can be fully adapted with a range of floats, wedges and raises to meet your clinical requirements. All Made to Order Footwear can have a maximum 3 upper adaptations per last.

All footwear contained in this catalogue is distributed by Blatchford in accordance with supplier agreements. Mendivil is a trade name of Biomedical footwear who have a designated (UKRP) United Kingdom Responsible Person for product compliance.

Contents

| | Page |
|---|-------|
| Mendivil Paediatric Footwear | 1 |
| Stocked Footwear | 3 |
| KSA Lazarus Boots and Bars | 5 |
| Made to Order Footwear | 7 |
| Specification Guide/Adaptations and Insoles | 7-8 |
| Colour Options | 9 |
| Ordering Information | 10 |
| Styles | 11-24 |
| Measurement Charts | 25-26 |
| Orthoses Accommodating Footwear | 27-28 |

The three Stocked Footwear styles on these pages are available for 24 hour delivery in the option colours and fastening shown.



Slickers Boot II



Sizes 27-45
Last Medium and Wide
Fastening Lace
Height UK
Stiffeners Standard
Colour Natural Leather Black

Part No. **Medium** 78001BK*M
Wide 78002BK*W * = Insert size required



Catapult Boot II



Sizes 27-45
Last Medium and Wide
Fastening Lace
Height UK
Stiffeners Standard
Colour Nubuck Black

Part No. **Medium** 78021BK*M
Wide 78022BK*W * = Insert size required



Sprinter Boot



Sizes 27-45
Last Medium and Wide
Fastening Lace
Height UK
Stiffeners Standard
Colour Nubuck Black

Part No. **Medium** 78501BK*M
Wide 78502BK*W * = Insert size required



Measurement Chart

| Size | Medium | | | | | Wide | | | | |
|------|--------|-------------|--------------|-------------|-----------|-------------|--------------|-------------|-----------|--|
| | Length | Joint Width | Instep Circ. | Ankle Circ. | Toe Depth | Joint Width | Instep Circ. | Ankle Circ. | Toe Depth | |
| 27 | 179 | 65 | 174 | 166 | 21 | 68 | 182 | 172 | 22 | |
| 28 | 186 | 67 | 178 | 171 | 22 | 70 | 186 | 177 | 22 | |
| 29 | 193 | 68 | 182 | 176 | 22 | 71 | 190 | 182 | 23 | |
| 30 | 199 | 69 | 186 | 180 | 23 | 72 | 194 | 188 | 23 | |
| 31 | 206 | 70 | 190 | 185 | 23 | 73 | 198 | 193 | 24 | |
| 32 | 213 | 72 | 194 | 190 | 24 | 75 | 202 | 198 | 25 | |
| 33 | 219 | 74 | 198 | 195 | 25 | 77 | 206 | 203 | 25 | |
| 34 | 226 | 75 | 202 | 200 | 25 | 78 | 210 | 208 | 26 | |
| 35 | 232 | 77 | 206 | 205 | 26 | 80 | 214 | 218 | 26 | |
| 36 | 239 | 79 | 210 | 210 | 26 | 82 | 218 | 224 | 27 | |
| 37 | 246 | 80 | 214 | 215 | 27 | 83 | 222 | 229 | 28 | |
| 38 | 252 | 82 | 218 | 220 | 28 | 85 | 226 | 235 | 28 | |
| 39 | 259 | 83 | 222 | 225 | 28 | 86 | 230 | 240 | 29 | |
| 40 | 266 | 85 | 226 | 230 | 29 | 88 | 234 | 246 | 29 | |
| 41 | 272 | 87 | 230 | 238 | 29 | 90 | 238 | 253 | 30 | |
| 42 | 279 | 88 | 234 | 246 | 30 | 91 | 242 | 259 | 31 | |
| 43 | 286 | 90 | 238 | 254 | 31 | 93 | 246 | 266 | 31 | |
| 44 | 292 | 92 | 242 | 252 | 31 | 95 | 250 | 272 | 32 | |
| 45 | 299 | 93 | 246 | 270 | 32 | 96 | 254 | 278 | 32 | |

Measurements for all sections are shown in mm.

KSA Lazarus Boots and Bars



| Plantar Flexion Stop Size | Product Code |
|---------------------------|--------------|
| 0000 | LP.PFS.0000 |
| 000 | LP.PFS.000 |
| 00 | LP.PFS.00 |
| 0 | LP.PFS.0 |
| 1 | LP.PFS.1 |
| 2 | LP.PFS.2 |
| 3 | LP.PFS.3 |
| 4 | LP.PFS.4 |
| 5 | LP.PFS.5 |
| 6 | LP.PFS.6 |
| 7 | LP.PFS.7 |
| 8 | LP.PFS.8 |

Lazarus Boot 2.0 - Grey

| Size | Product Code Left | Product Code Right | Product Code Pair |
|------|-------------------|--------------------|-------------------|
| 6X0 | LP.GR.S.000000.L | LP.GR.S.000000.R | LP.GR.S.000000 |
| 5X0 | LP.GR.S.00000.L | LP.GR.S.00000.R | LP.GR.S.00000 |
| 0000 | LP.GR.S.0000.L | LP.GR.S.0000.R | LP.GR.S.0000 |
| 000 | LP.GR.S.000.L | LP.GR.S.000.R | LP.GR.S.000 |
| 00 | LP.GR.S.00.L | LP.GR.S.00.R | LP.GR.S.00 |
| 0 | LP.GR.S.0.L | LP.GR.S.0.R | LP.GR.S.0 |
| 1 | LP.GR.S.1.L | LP.GR.S.1.R | LP.GR.S.1 |
| 2 | LP.GR.S.2.L | LP.GR.S.2.R | LP.GR.S.2 |
| 3 | LP.GR.S.3.L | LP.GR.S.3.R | LP.GR.S.3 |
| 4 | LP.GR.S.4.L | LP.GR.S.4.R | LP.GR.S.4 |
| 5 | LP.GR.S.5.L | LP.GR.S.5.R | LP.GR.S.5 |
| 6 | LP.GR.S.6.L | LP.GR.S.6.R | LP.GR.S.6 |
| 7 | LP.GR.S.7.L | LP.GR.S.7.R | LP.GR.S.7 |
| 8 | LP.GR.S.8.L | LP.GR.S.8.R | LP.GR.S.8 |
| 9 | LP.GR.S.9.L | LP.GR.S.9.R | LP.GR.S.9 |
| 10 | LP.GR.S.10.L | LP.GR.S.10.R | LP.GR.S.10 |
| 11 | LP.GR.S.11.L | LP.GR.S.11.R | LP.GR.S.11 |
| 12 | LP.GR.S.12.L | LP.GR.S.12.R | LP.GR.S.12 |

Lazarus Boot 2.0 - Royal Blue

| Size | Product Code Left | Product Code Right | Product Code Pair |
|------|-------------------|--------------------|-------------------|
| 6X0 | LP.RB.S.000000.L | LP.RB.S.000000.R | LP.RB.S.000000 |
| 5X0 | LP.RB.S.00000.L | LP.RB.S.00000.R | LP.RB.S.00000 |
| 0000 | LP.RB.S.0000.L | LP.RB.S.0000.R | LP.RB.S.0000 |
| 000 | LP.RB.S.000.L | LP.RB.S.000.R | LP.RB.S.000 |
| 00 | LP.RB.S.00.L | LP.RB.S.00.R | LP.RB.S.00 |
| 0 | LP.RB.S.0.L | LP.RB.S.0.R | LP.RB.S.0 |
| 1 | LP.RB.S.1.L | LP.RB.S.1.R | LP.RB.S.1 |
| 2 | LP.RB.S.2.L | LP.RB.S.2.R | LP.RB.S.2 |
| 3 | LP.RB.S.3.L | LP.RB.S.3.R | LP.RB.S.3 |
| 4 | LP.RB.S.4.L | LP.RB.S.4.R | LP.RB.S.4 |
| 5 | LP.RB.S.5.L | LP.RB.S.5.R | LP.RB.S.5 |
| 6 | LP.RB.S.6.L | LP.RB.S.6.R | LP.RB.S.6 |
| 7 | LP.RB.S.7.L | LP.RB.S.7.R | LP.RB.S.7 |
| 8 | LP.RB.S.8.L | LP.RB.S.8.R | LP.RB.S.8 |
| 9 | LP.RB.S.9.L | LP.RB.S.9.R | LP.RB.S.9 |
| 10 | LP.RB.S.10.L | LP.RB.S.10.R | LP.RB.S.10 |
| 11 | LP.RB.S.11.L | LP.RB.S.11.R | LP.RB.S.11 |
| 12 | LP.RB.S.12.L | LP.RB.S.12.R | LP.RB.S.12 |

Lazarus Boot 2.0 - Pink

| Size | Product Code Left | Product Code Right | Product Code Pair |
|------|-------------------|--------------------|-------------------|
| 6X0 | LP.PI.S.000000.L | LP.PI.S.000000.R | LP.PI.S.000000 |
| 5X0 | LP.PI.S.00000.L | LP.PI.S.00000.R | LP.PI.S.00000 |
| 0000 | LP.PI.S.0000.L | LP.PI.S.0000.R | LP.PI.S.0000 |
| 000 | LP.PI.S.000.L | LP.PI.S.000.R | LP.PI.S.000 |
| 00 | LP.PI.S.00.L | LP.PI.S.00.R | LP.PI.S.00 |
| 0 | LP.PI.S.0.L | LP.PI.S.0.R | LP.PI.S.0 |
| 1 | LP.PI.S.1.L | LP.PI.S.1.R | LP.PI.S.1 |
| 2 | LP.PI.S.2.L | LP.PI.S.2.R | LP.PI.S.2 |
| 3 | LP.PI.S.3.L | LP.PI.S.3.R | LP.PI.S.3 |
| 4 | LP.PI.S.4.L | LP.PI.S.4.R | LP.PI.S.4 |
| 5 | LP.PI.S.5.L | LP.PI.S.5.R | LP.PI.S.5 |
| 6 | LP.PI.S.6.L | LP.PI.S.6.R | LP.PI.S.6 |
| 7 | LP.PI.S.7.L | LP.PI.S.7.R | LP.PI.S.7 |
| 8 | LP.PI.S.8.L | LP.PI.S.8.R | LP.PI.S.8 |
| 9 | LP.PI.S.9.L | LP.PI.S.9.R | LP.PI.S.9 |
| 10 | LP.PI.S.10.L | LP.PI.S.10.R | LP.PI.S.10 |
| 11 | LP.PI.S.11.L | LP.PI.S.11.R | LP.PI.S.11 |
| 12 | LP.PI.S.12.L | LP.PI.S.12.R | LP.PI.S.12 |



Sherling Pad - Universal One Size

| Product Code | Description |
|--------------|--|
| LP.ZPT.GR | Grey Shearling Padded Tongue Pack 10 |
| LP.ZPT.PI | Pink Shearling Padded Tongue Pack 10 |
| LP.ZPT.BU | Burgundy Shearling Padded Tongue Pack 10 |
| LP.ZPT.RB | Royal Blue Shearling Padded Tongue Pack 10 |
| LP.ZPT.NB | Navy Blue Shearling Padded Tongue Pack 10 |

- All bars are fully assembled
- Simple quick release mechanism
- Set the width of the bar and the angle to your patient preference
- Fasten the boots to the bar with the click system.



Product Code: LPSB.0.S
Static Bar Small - 0 Degree



Product Code: LPSB.0.L
Static Bar Large - 0 Degree



Product Code: LPSB.15.S
Static Bar Small - 15 Degree



Product Code: LPSB.15.L
Static Bar Large - 15 Degree



Product Code: LPAB.0.S
Dynamic Bar Small



Product Code: LPAB.0.L
Dynamic Bar Large



| Size | Length | Width |
|--------|-----------|---------|
| 000000 | 5.5-6.5 | 2.5-4.0 |
| 00000 | 5.8-6.9 | 3.0-4.5 |
| 0000 | 6.0-7.3 | 3.3-4.8 |
| 000 | 6.3-7.8 | 3.0-5.5 |
| 00 | 6.8-8.3 | 3.0-5.5 |
| 0 | 7.5-9.0 | 3.0-6.0 |
| 1 | 8.0-9.5 | 3.5-6.0 |
| 2 | 9.0-10.5 | 3.8-6.5 |
| 3 | 10.0-11.5 | 4.5-7.0 |
| 4 | 11.0-12.5 | 4.7-7.0 |


| Size | Length | Width |
|-----------|-------------|---------|
| 5 | 12.0-13.5 | 4.8-7.0 |
| 6 | 13.0-14.5 | 5.0-7.5 |
| 7 | 14.0-15.5 | 5.5-8.0 |
| 8 | 15.0-16.5 | 6.7-8.2 |
| 9 | 16.0-17.5 | 7.0-8.5 |
| 10 | 17.0-18.5 | 7.4-8.9 |
| 11 | 18.0-19.5 | 7.7-9.2 |
| 12 | 19.0-20.5 | 8.0-9.5 |
| Small Bar | 145mm-212mm | |
| Large Bar | 217mm-355mm | |



Made to Order Specification Guide

Please clearly state the following options when ordering your chosen style. Not all options are applicable to every style, please check the specification provided alongside each boot in the catalogue pages.

- Choose from either a **lace** or **reverse hook and loop** fastening (hook and loop incurs an additional charge).
- Specify the correct height:



- 5 **standard**, 3 **anti-varus** and 3 **orthoses accommodating last** options are available. See page 25 for measurement charts. (Split sizes are available but will incur an additional charge).
- 3 **stiffener** options are available:




- Last additions: a maximum of 3 additions per last from the 7 illustrated** are available with modifications in multiples of 3mm, allowing you to choose any combination as illustrated below. Total modifications must not exceed 18mm per pair, each modification is charged as one unit with the code LADP-PAEDS.

If the limits are exceeded then a fully drafted Made To Measure order is required. If the footwear is to be size 32 or smaller, a cast is also required.

All footwear is supplied complete fittings are not available.

Last adaptations are not available in Accommodating Orthoses footwear.



- 5 **sole options** are available in the size ranges shown below:



- Various **colour and material options** can be specified for each style, as illustrated on page 7.

- All 5 **standard last width MTO footwear** is supplied with the following internal inlay.



Adaptations and Insoles


| | |
|-------------|--|
| A17-PAEDS | High boot uppers additional 25mm to footwear |
| B1 | Outside compensating raise up to 25mm with a rand |
| B2 | Each subsequent 25mm or part thereof to B1 |
| B3 | Outside compensating raise up to 25mm without a rand |
| B4 | Each subsequent 25mm or part thereof to B3 |
| B5 | Raising the heel in a single item of footwear up to 25mm |
| B7 | Addition of a heel seat or valgus/varus wedge fitted to shoe (each) |
| B8-PAEDS | Stirrup, heel or caliper retaining strap (each) |
| B9 | Any additional padding in a single item of footwear |
| B15-PAEDS | Addition of orthosis stiffener |
| XSTRG-PAEDS | Extra Strength Stability Stiffener |
| B16 | Addition of a metatarsal bar fitted to the sole (each) |
| B17 | Inclusion of a rocker sole (each) |
| B18 | Inclusion of wedge, fitted into/onto the sole (each) |
| B19 | Inclusion of a floated sole in a single item of footwear |
| B20 | Inclusion of a sole strip or plate in composite material (each) |
| B21 | Inclusion of a sole strip or plate in non-composite material (each) |
| B22 | Addition of a rubber tip to the sole (each) |
| B23 | Addition of a tip, with flange, to protect the welt (each) |
| B24 | Addition of an external protective toe cap (each) |
| B25 | Full/half-elongated, oblique, floated heel or heel wedge or torque heel in a single item of footwear |
| B30 | Socket to one/both sides of the heel, without front/backstop to single item of footwear |
| B31 | Socket to one/both sides of the heel, with front/backstop to single item of footwear |
| B32 | Inclusion of a toe raising socket to single item of footwear |
| B34-PAEDS | Addition of a lined and padded double-decker T or Y strap (each) |
| B35 | Rebuilding the heel to allow socket fixation in a single item of footwear |
| B40a | 3/4 or full length non-moulded inlay, non-shock absorb material e.g. Leather/Neoprene/Plastazote® (each) |
| B40b | 3/4 or full length non-moulded inlay, shock absorb material e.g. Poron®/PPT® (each) |
| B40c | 3/4 or full length non-moulded inlay in silicone type material (each) |
| B41 | Metatarsal/valgus or heel pad to inlay (each) |
| B42 | Metatarsal/valgus or heel pad fitted direct to shoe (each) |
| B43 | Inclusion of a heel elevator in cork, HDP or EVA covered (each) |
| B50 | Filler to compensate for toe amputation(s) (each) |
| B51 | Filler to compensate for amputation of forepart of foot (each) |
| LADP-PAEDS | Last supplement per adjustment (each) |
| NSB030 | Flat carbon fibre insole (each) |
| PEADHL | Paediatric Footwear H&L Fastening Option |
| SURMENDZR | Paediatric Footwear Split Size Charge up to Size 38 |
| SURMENDS | Paediatric Footwear Split Size Charge Over Size 39 |
| WLN-PAEDS | Addition of lambswool lining in a single item of footwear (each) |

(The dark blue lines on the diagram to the right illustrate the area lined.)



Colour Options

Not all colour options are applicable to every style, please check the specification provided alongside each boot in the catalogue pages.




Leather - Natural
Colour - Black



Leather - Natural
Colour - Navy




Leather - Nubuck
Colour - Tartufo



Leather - Box
Colour - White



Leather - Nubuck
Colour - Black




Leather - Nubuck
Colour - Chine




Leather - Val
Colour - Linen



Leather - Nubuck
Colour - Pink



Leather - Patent
Colour - Black



Leather - Natural
Colour - Royal




Leather - Nubuck
Colour - Light Brown



Leather - Bolonia
Colour - Pink



Leather - Calgary
Colour - Black




Leather - Nubuck
Colour - Navy




Leather - Val
Colour - Soan




Leather - Nubuck
Colour - Fuchsia



Leather - Nubuck
Colour - Grey



Leather - Antic
Colour - Jean



Leather - Val
Colour - Saddale



Leather - Bolonia
Colour - Fuchsia



Leather - Bolonia
Colour - Lilac

Ordering Information

1. Ordering and Enquiries:

Ordering Paediatric footwear from Blatchford is simple and easy. You can contact our Customer Services directly:

Telephone: +44 (0) 114 2637900

Email: orthoticscs@blatchford.co.uk

Facsimile: +44 (0) 114 2637901

Post: Blatchford Limited, 11 Atlas Way, Atlas North, Sheffield, S4 7QQ.

2. Delivery Promise:

Our goal is to provide fast and efficient delivery. We aim to despatch all orders for Paediatric Made to Order footwear from receipt in 14 working days. This does not apply where specific delivery schedules have been previously agreed.

3. Catalogues/Product Information:

If you require additional catalogues or product information please contact our customer service advisors who will be happy to help with your enquiry.

4. Single Patient Use:

All footwear featured in this catalogue is prescribed/distributed by Blatchford for single patient use only and should be fitted by a qualified healthcare professional. The reuse of any footwear can potentially cause harm to the user and will invalidate the UKCA/CE or any warranty agreement. Please refer to the IFU for more information.

5. Shortages/Damaged Goods:

If any items are received damaged or there are shortages in the delivery this must be notified in writing to customer services within 5 working days of receipt of goods.

6. Returns:

In the unlikely event of any footwear being supplied faulty or having been incorrectly supplied, please contact customer services so that the correct goods are despatched immediately. We are happy to accept Paediatric Stock Footwear that is returned for credit if they are not suitable for the patient.

For items to be accepted for credit the following conditions apply:

- The footwear must be unworn and returned in less than 30 days from the date of the delivery
- The product must not have been adapted in any way
- Split sized or made to order footwear is non returnable for credit
- Split sized footwear direct from stock is returnable for credit
- Are returned in the original packaging
- Are in a re-saleable condition

Blatchford reserves the right to reject items returned for credit if the above conditions are not met. Blatchford is unable to accept return of items that are 'made to order' for credit, unless the goods are faulty or it can be demonstrated they have been damaged in transit. Products and prices are subject to variation at any time and are subject to availability. E & OE.

Need boots for school? No problem!

All of the styles shown in this catalogue are available in a range of black leathers, making them perfect for school days.

Shikra




Last Options 5 standard widths and 3 orthoses widths
Sizes 16-44 (subject to sole option)
Anti-Varus Last 3 widths - sizes 16-33 (subject to sole option)
Fastening Lace or Reverse Pull Hook and Loop (with medial zip)
Height UK
Stiffeners Standard, Orthoses Accommodation

PU Sole (Sizes 16-45)
 Black White Beige Grey Pink

Rubber Kick (Sizes 20-45)
 Black Light Beige




Not suitable with extra strength stability stiffener.



Leather - Bolonia Colour - Pink

Leather - Natural Colour - Royal



Leather - Bolonia Colour - Lilac

Leather - Nubuck Colour - Grey



Leather - Nubuck Colour - Chine

Leather - Nubuck Colour - Tartufo


Owl



Last Options 5 standard widths and 3 orthoses widths
Sizes 16-45
Anti-Varus Last 3 widths - sizes 16-33
Fastening Lace or Reverse Pull Hook and Loop (with medial zip)
Height UK
Stiffeners Standard, Orthoses Accommodation

PU Sole (Sizes 16-45)
 Black White Beige Grey Pink

Not suitable with extra strength stability stiffener.




Leather - Nubuck Colour - Fuchsia

Leather - Nubuck Colour - Black



Leather - Nubuck Colour - Chine

Leather - Nubuck Colour - Tartufo



Leather - Natural Colour - Black

Leather - Natural Colour - Royal

Raptor



Last Options 5 standard widths and 3 orthoses widths
Sizes 16-45 (subject to sole option)
Anti-Varus Last 3 widths - sizes 16-33 (subject to sole option)
Fastening Lace or Reverse Pull Hook and Loop
Height UK
Stiffeners Standard, Orthoses Accommodation and Extra Strength Stability


PU Sole (Sizes 16-45)
 Black White Beige Grey Pink

Rubber Montana (Sizes 18-45)
 Black Beige



Leather - Bolonia Colour - Lilac

Leather - Calgary Colour - Black



Leather - Nubuck Colour - Tartufo

Leather - Natural Colour - Black



Leather - Natural Colour - Royal

Eagle




Last Options 5 standard widths and 3 orthoses widths
Sizes 16-44 (subject to sole option)
Anti-Varus Last 3 widths - sizes 16-33 (subject to sole option)
Fastening Lace or Reverse Pull Hook and Loop
Height UK
Stiffeners Standard, Orthoses Accommodation and Extra Strength Stability

PU Sole (Sizes 16-45)
 Black White Beige Grey Pink


Rubber Montana (Sizes 18-45)
 Black Beige

Rubber Kick (Sizes 20-45)
 Black Light Beige




Leather - Nubuck Colour - Pink

Leather - Calgary Colour - Black



Leather - Nubuck Colour - Grey

Leather - Bolonia Colour - Fuchsia



Leather - Nubuck Colour - Fuchsia

Kite



Last Options 5 standard widths and 3 orthoses widths
Sizes 16-44 (subject to sole option)
Anti-Varus Last 3 widths - sizes 16-33 (subject to sole option)
Fastening Lace or Reverse Pull Hook and Loop
Height UK
Stiffeners Standard, Orthoses Accommodation and Extra Strength Stability



Osprey




Last Options 5 standard widths and 3 orthoses widths
Sizes 16-44 (subject to sole option)
Anti-Varus Last 3 widths - sizes 16-33 (subject to sole option)
Fastening Lace or Reverse Pull Hook and Loop
Height UK or Continental
Stiffeners Standard, Orthoses Accommodation and Extra Strength Stability



Condor



Last Options 5 standard widths and 3 orthoses widths
Sizes 18-45
Anti-Varus Last 3 widths - sizes 18-33
Fastening Lace or Reverse Pull Hook and Loop (with medial zip)
Height Continental
Stiffeners Standard, Orthoses Accommodation



Harrier



Last Options 5 standard widths and 3 orthoses widths
Sizes 16-45 (subject to sole option)
Anti-Varus Last 3 widths - sizes 16-33 (subject to sole option)
Fastening Lace or Reverse Pull Hook and Loop
Height UK or Continental
Stiffeners Standard, Orthoses Accommodation and Extra Strength Stability



Bateleur



Last Options 5 standard widths and 3 orthoses widths
Sizes 16-45 (subject to sole option)
Anti-Varus Last 3 widths - sizes 16-33 (subject to sole option)
Fastening Lace or Reverse Pull Hook and Loop
Height UK
Stiffeners Standard, Orthoses Accommodation and Extra Strength Stability




Gryfalcon



Last Options 5 standard widths and 3 orthoses widths
Sizes 16-45 (subject to sole option)
Anti-Varus Last 3 widths - sizes 16-33 (subject to sole option)
Fastening Lace or Reverse Pull Hook and Loop
Height UK
Stiffeners Standard, Orthoses Accommodation and Extra Strength Stability




Only available in colour options below.



Kestrel



Last Options 5 standard widths and 3 orthoses widths
Sizes 16-45 (subject to sole option)
Anti-Varus Last 3 widths - sizes 16-33 (subject to sole option)
Fastening Lace or Reverse Pull Hook and Loop (with medial zip)
Height UK or Continental
Stiffeners Standard, Orthoses Accommodation



Kookaburra



Last Options 5 standard widths and 3 orthoses widths
Sizes 16-45 (subject to sole option)
Anti-Varus Last 3 widths - sizes 16-33 (subject to sole option)
Fastening Lace or Reverse Pull Hook and Loop
Height UK or Continental
Stiffeners Standard, Orthoses Accommodation and Extra Strength Stability



Only available in colour options below.



Merlin



Last Options 5 standard widths and 3 orthoses widths
Sizes 16-45
Anti-Varus Last 3 widths - sizes 16-33
Fastening Lace or Reverse Pull Hook and Loop
Height UK
Stiffeners Standard, Orthoses Accommodation and Extra Strength Stability



Only available in colour options below.



Peregrine



Last Options 5 standard widths and 3 orthoses widths
Sizes 16-45
Anti-Varus Last 3 widths - sizes 16-33
Fastening Lace or Reverse Pull Hook and Loop
Height UK or Continental
Stiffeners Standard, Orthoses Accommodation and Extra Strength Stability




Only available in colour options below.



Besra



Last Options 5 standard widths and 3 orthoses widths
Sizes 16-45
Anti-Varus Last 3 widths - sizes 16-33
Fastening Lace or Reverse Pull Hook and Loop
Height UK
Stiffeners Standard, Orthoses Accommodation and Extra Strength Stability



Bazza



Last Options 5 standard widths and 3 orthoses widths
Sizes 16-45
Anti-Varus Last 3 widths - sizes 16-33
Fastening Lace or Reverse Pull Hook and Loop
Height UK or Continental
Stiffeners Standard, Orthoses Accommodation and Extra Strength Stability



Hobby



Last Options 5 standard widths and 3 orthoses widths
Sizes 16-45
Anti-Varus Last 3 widths - sizes 16-33
Fastening Lace or Reverse Pull Hook and Loop
Height UK or Continental
Stiffeners Standard, Orthoses Accommodation and Extra Strength Stability



Goshawk

Black rubber heel and toe protector is standard on design.



Last Options 5 standard widths and 3 orthoses widths
Sizes 16-45 (subject to sole option)
Anti-Varus Last 3 widths - sizes 16-33 (subject to sole option)
Fastening Lace or Reverse Pull Hook and Loop
Height UK or Continental
Stiffeners Standard, Orthoses Accommodation and Extra Strength Stability



Only available in colour options below.



Buzzard



Last Options 5 standard widths and 3 orthoses widths
Sizes 16-45 (subject to sole option)
Anti-Varus Last 3 widths - sizes 16-33 (subject to sole option)
Fastening Lace or Reverse Pull Hook and Loop
Height UK
Stiffeners Standard, Orthoses Accommodation and Extra Strength Stability



Jay



Last Options 5 standard widths and 3 orthoses widths
Sizes 16-45 (subject to sole option)
Anti-Varus Last 3 widths - sizes 16-33
Fastening Lace only with medial zip
Height UK or Continental
Stiffeners Standard, Orthoses Accommodation



Only available in colour options below.



Eider



Last Options 5 standard widths and 3 orthoses widths
Sizes 16-45 (subject to sole option)
Anti-Varus Last 3 widths - sizes 16-33 (subject to sole option)
Fastening Lace or Reverse Pull Hook and Loop
Height UK
Stiffeners Standard, Orthoses Accommodation and Extra Strength Stability



Raven



Last Options 5 standard widths and 3 orthoses widths
Sizes 16-45 (subject to sole option)
Anti-Varus Last 3 widths - sizes 16-33 (subject to sole option)
Fastening Lace or Reverse Pull Hook and Loop
Height UK
Stiffeners Standard, Orthoses Accommodation and Extra Strength Stability



Cub Starter Boot



Last Options 5 standard widths
Sizes 16-26
Anti-Varus Last 3 widths - sizes 16-26
Fastening Lace or Reverse Pull Hook and Loop
Height UK
Stiffeners Standard, Orthoses Accommodation and Extra Strength Stability



Slickers Boot II



Last Options 5 standard widths
Sizes 18-45
Anti-Varus Last 3 widths - sizes 18-33
Fastening Lace or Reverse Pull Hook and Loop
Height UK
Stiffeners Standard, Orthoses Accommodation and Extra Strength Stability



Sprinter Boot



Last Options 5 standard widths
Sizes 16-45
Anti-Varus Last 3 widths - sizes 16-33
Fastening Lace or Reverse Pull Hook and Loop
Height UK
Stiffeners Standard, Orthoses Accommodation and Extra Strength Stability



Catapult Boot II



Last Options 5 standard widths
Sizes 18-45
Anti-Varus Last 3 widths - sizes 18-33
Fastening Lace or Reverse Pull Hook and Loop
Height UK
Stiffeners Standard, Orthoses Accommodation and Extra Strength Stability



Saker Sandal



Last Options 5 standard widths and 3 orthoses widths
Sizes 16-45
Fastening Reverse Pull Hook and Loop
Height Continental
Stiffeners Standard, Orthoses Accommodation



Not suitable with extra strength stability stiffener.
Only available in colour options below.




Measurement Charts

Fitting Principles

The measurement charts lay out a step by step procedure which can help the clinical assessment and measurement process achieve a right first time fit.

It is important to note that when the foot weight bears, the foot posture changes and can have an effect on the shoe size required. Please consider this to ensure the desired toe room is added to the overall length.



1. Foot length
2. Joint width
3. Instep circumference
4. Ankle circumference

Standard Last

| Boot | | Narrow | | | | Medium | | | | Wide | | | | X Wide | | | | XX Wide | | | |
|------|--------|-------------|--------------|-------------|-----------|-------------|--------------|-------------|-----------|-------------|--------------|-------------|-----------|-------------|--------------|-------------|-----------|-------------|--------------|-------------|-----------|
| Size | Length | Joint Width | Instep Circ. | Ankle Circ. | Toe Depth | Joint Width | Instep Circ. | Ankle Circ. | Toe Depth | Joint Width | Instep Circ. | Ankle Circ. | Toe Depth | Joint Width | Instep Circ. | Ankle Circ. | Toe Depth | Joint Width | Instep Circ. | Ankle Circ. | Toe Depth |
| 16 | 106 | 46 | 122 | 104 | 14 | 49 | 130 | 106 | 14 | 52 | 138 | 119 | 15 | 55 | 146 | 120 | 15 | 58 | 154 | 122 | 16 |
| 17 | 113 | 47 | 126 | 108.5 | 15 | 50 | 134 | 111.5 | 15 | 53 | 142 | 124 | 16 | 56 | 150 | 124.5 | 16 | 59 | 158 | 127.5 | 17 |
| 18 | 119 | 48 | 130 | 113 | 15 | 51 | 138 | 117 | 16 | 54 | 146 | 129 | 16 | 57 | 154 | 129 | 17 | 60 | 162 | 133 | 17 |
| 19 | 126 | 50 | 134 | 117.5 | 16 | 53 | 142 | 122 | 16 | 56 | 150 | 134 | 17 | 59 | 158 | 133.5 | 17 | 62 | 166 | 138 | 18 |
| 20 | 133 | 51 | 138 | 122 | 16 | 54 | 146 | 128 | 17 | 57 | 154 | 138 | 17 | 60 | 162 | 138 | 18 | 63 | 170 | 144 | 18 |
| 21 | 139 | 53 | 142 | 128 | 17 | 56 | 150 | 133.5 | 17 | 59 | 158 | 143 | 18 | 62 | 166 | 143.5 | 18 | 65 | 174 | 149 | 19 |
| 22 | 146 | 54 | 146 | 134 | 18 | 57 | 154 | 139 | 18 | 60 | 162 | 148 | 19 | 63 | 170 | 149 | 19 | 66 | 178 | 154 | 20 |
| 23 | 153 | 55 | 150 | 140 | 18 | 58 | 158 | 144.5 | 19 | 61 | 166 | 153 | 19 | 64 | 174 | 154.5 | 20 | 67 | 182 | 159 | 20 |
| 24 | 159 | 56 | 154 | 146 | 19 | 59 | 162 | 150 | 19 | 62 | 170 | 158 | 20 | 65 | 178 | 160 | 20 | 68 | 186 | 164 | 21 |
| 25 | 166 | 58 | 158 | 150 | 19 | 61 | 166 | 156 | 20 | 64 | 174 | 162 | 20 | 67 | 182 | 166 | 21 | 70 | 190 | 168 | 21 |
| 26 | 173 | 60 | 162 | 154 | 20 | 63 | 170 | 161 | 20 | 66 | 178 | 167 | 21 | 69 | 186 | 171 | 21 | 72 | 194 | 173 | 22 |
| 27 | 179 | 62 | 166 | 158 | 21 | 65 | 174 | 166 | 21 | 68 | 182 | 172 | 22 | 71 | 190 | 176 | 22 | 74 | 198 | 178 | 23 |
| 28 | 186 | 64 | 170 | 162 | 21 | 67 | 178 | 171 | 22 | 70 | 186 | 177 | 22 | 73 | 194 | 181 | 23 | 76 | 202 | 183 | 23 |
| 29 | 193 | 65 | 174 | 166 | 22 | 68 | 182 | 176 | 22 | 71 | 190 | 182 | 23 | 74 | 198 | 186 | 23 | 77 | 206 | 188 | 24 |
| 30 | 199 | 66 | 178 | 170 | 22 | 69 | 186 | 180 | 23 | 72 | 194 | 188 | 23 | 75 | 202 | 190 | 24 | 78 | 210 | 193 | 24 |
| 31 | 206 | 67 | 182 | 175 | 23 | 70 | 190 | 185 | 23 | 73 | 198 | 193 | 24 | 76 | 206 | 196 | 24 | 79 | 214 | 198 | 25 |
| 32 | 213 | 69 | 186 | 180 | 24 | 72 | 194 | 190 | 24 | 75 | 202 | 198 | 25 | 78 | 210 | 202 | 25 | 81 | 218 | 205 | 26 |
| 33 | 219 | 71 | 190 | 185 | 24 | 74 | 198 | 195 | 25 | 77 | 206 | 203 | 25 | 80 | 214 | 208 | 26 | 83 | 222 | 212 | 26 |
| 34 | 226 | 72 | 194 | 190 | 25 | 75 | 202 | 200 | 25 | 78 | 210 | 208 | 26 | 81 | 218 | 214 | 26 | 84 | 226 | 219 | 27 |
| 35 | 232 | 74 | 198 | 196 | 25 | 77 | 206 | 205 | 26 | 80 | 214 | 218 | 26 | 83 | 222 | 220 | 27 | 86 | 230 | 226 | 27 |
| 36 | 239 | 76 | 202 | 201.5 | 26 | 79 | 210 | 210 | 26 | 82 | 218 | 223.5 | 27 | 85 | 226 | 226.5 | 27 | 88 | 234 | 231.5 | 28 |
| 37 | 246 | 77 | 206 | 207 | 27 | 80 | 214 | 215 | 27 | 83 | 222 | 229 | 28 | 86 | 230 | 232 | 28 | 89 | 238 | 237 | 29 |
| 38 | 252 | 79 | 210 | 212.5 | 27 | 82 | 218 | 220 | 28 | 85 | 226 | 234.5 | 28 | 88 | 234 | 238 | 29 | 91 | 242 | 242.5 | 29 |
| 39 | 259 | 80 | 214 | 218.5 | 28 | 83 | 222 | 225 | 28 | 86 | 230 | 240 | 29 | 89 | 238 | 244 | 29 | 92 | 246 | 248 | 30 |
| 40 | 266 | 82 | 218 | 224 | 28 | 85 | 226 | 230 | 29 | 88 | 234 | 246 | 29 | 91 | 242 | 252 | 30 | 94 | 250 | 253.5 | 30 |
| 41 | 272 | 84 | 222 | 233 | 29 | 87 | 230 | 238 | 29 | 90 | 238 | 252.5 | 30 | 93 | 246 | 257.5 | 30 | 96 | 254 | 260 | 31 |
| 42 | 279 | 85 | 226 | 242 | 30 | 88 | 234 | 246 | 30 | 91 | 242 | 259 | 31 | 94 | 250 | 263 | 31 | 97 | 258 | 268.5 | 32 |
| 43 | 286 | 87 | 230 | 252 | 30 | 90 | 238 | 254 | 31 | 93 | 246 | 265.5 | 31 | 96 | 254 | 268.5 | 32 | 99 | 262 | 277 | 32 |
| 44 | 292 | 89 | 234 | 260 | 31 | 92 | 242 | 252 | 31 | 95 | 250 | 272 | 32 | 98 | 258 | 274 | 32 | 101 | 266 | 285.5 | 33 |
| 45 | 299 | 90 | 238 | 268 | 31 | 93 | 246 | 270 | 32 | 96 | 254 | 278 | 32 | 99 | 262 | 280 | 33 | 102 | 270 | 294 | 33 |

Measurements for all sections are shown in mm.



Anti-Varus Last

| Boot | | Medium | | | Wide | | | X Wide | | |
|------|--------|-------------|--------------|-------------|-------------|--------------|-------------|-------------|--------------|-------------|
| Size | Length | Joint Width | Instep Circ. | Ankle Circ. | Joint Width | Instep Circ. | Ankle Circ. | Joint Width | Instep Circ. | Ankle Circ. |
| 16 | 106 | 49 | 130 | 106 | 52 | 138 | 119 | 55 | 146 | 120 |
| 17 | 113 | 50 | 134 | 111.5 | 53 | 142 | 124 | 56 | 150 | 124.5 |
| 18 | 119 | 51 | 138 | 117 | 54 | 146 | 129 | 57 | 154 | 129 |
| 19 | 126 | 53 | 142 | 122 | 56 | 150 | 134 | 59 | 158 | 133.5 |
| 20 | 133 | 54 | 146 | 128 | 57 | 154 | 138 | 60 | 162 | 138 |
| 21 | 139 | 56 | 150 | 133.5 | 59 | 158 | 143 | 62 | 166 | 143.5 |
| 22 | 146 | 57 | 154 | 139 | 60 | 162 | 148 | 63 | 170 | 149 |
| 23 | 153 | 58 | 158 | 144.5 | 61 | 166 | 153 | 64 | 174 | 154.5 |
| 24 | 159 | 59 | 162 | 150 | 62 | 170 | 158 | 65 | 178 | 160 |
| 25 | 166 | 61 | 166 | 156 | 64 | 174 | 162 | 67 | 182 | 166 |
| 26 | 173 | 63 | 170 | 161 | 66 | 178 | 167 | 69 | 186 | 171 |
| 27 | 179 | 65 | 174 | 166 | 68 | 182 | 172 | 71 | 190 | 176 |
| 28 | 186 | 67 | 178 | 171 | 70 | 186 | 177 | 73 | 194 | 181 |
| 29 | 193 | 68 | 182 | 176 | 71 | 190 | 182 | 74 | 198 | 186 |
| 30 | 199 | 69 | 186 | 180 | 72 | 194 | 188 | 75 | 202 | 190 |
| 31 | 206 | 70 | 190 | 185 | 73 | 198 | 193 | 76 | 206 | 196 |
| 32 | 213 | 72 | 194 | 190 | 75 | 202 | 198 | 78 | 210 | 202 |
| 33 | 219 | 74 | 198 | 195 | 77 | 206 | 203 | 80 | 214 | 208 |

Measurements for all sections are shown in mm.

Orthoses Accommodating Footwear


The styles below feature an additional zip fastening for ease of access and donning.

Falcon



Last Options 3 orthoses widths
Sizes 20-44
Fastening Lace with Zip
Height UK
Stiffeners Standard, Orthoses Accommodation

No last adaptation available on Accommodating Footwear.


Rubber Montana 
 (Sizes 18-45) Black (with yellow stitch effect)

Rubber Kick 
 (Sizes 20-45) Black (with yellow stitch effect)

Not suitable with extra strength stability stiffener.



Leather - Calgary
Colour - Black Leather - Nubuck
Colour - Black



Leather - Natural
Colour - Black Leather - Natural
Colour - Royal



Leather - Nubuck
Colour - Brown

Hawk




Last Options 3 orthoses widths
Sizes 20-45
Fastening Lace or Reverse Pull Hook and Loop with Zip
Height UK
Stiffeners Standard, Orthoses Accommodation

No last adaptation available on Accommodating Footwear.

Sports 
 (Sizes 16-45) Black Beige Grey White

Not suitable with extra strength stability stiffener.



Leather - Box
Colour - White Leather - Natural
Colour - Black



Leather - Natural
Colour - Royal Leather - Bolonia
Colour - Red



Leather - Bolonia
Colour - Fuchsia Leather - Bolonia
Colour - Pink

Orthosis Accommodating Last

| Boot | | Orthoses Narrow Last | | | | Orthoses Wide Last | | | | Orthoses Extra Volume | | | |
|------|--------|----------------------|-------------|-------------|-----------|--------------------|-------------|-------------|-----------|-----------------------|-------------|-------------|-----------|
| Size | Length | Joint Width | Joint Circ. | Ankle Circ. | Long Heel | Joint Width | Joint Circ. | Ankle Circ. | Long Heel | Joint Width | Joint Circ. | Ankle Circ. | Long Heel |
| 20 | 133 | 57 | 170 | 161 | 225 | 60 | 178 | 168 | 235 | 60 | 184 | 168 | 241 |
| 21 | 139 | 59 | 174 | 167 | 233 | 62 | 182 | 174 | 243 | 62 | 188 | 174 | 249 |
| 22 | 146 | 60 | 178 | 173 | 241 | 63 | 186 | 180 | 251 | 63 | 192 | 180 | 257 |
| 23 | 153 | 61 | 182 | 179 | 249 | 64 | 190 | 186 | 259 | 64 | 196 | 186 | 265 |
| 24 | 159 | 62 | 186 | 185 | 257 | 65 | 194 | 192 | 267 | 65 | 200 | 192 | 273 |
| 25 | 166 | 64 | 190 | 191 | 265 | 67 | 198 | 198 | 275 | 67 | 204 | 198 | 281 |
| 26 | 173 | 66 | 194 | 197 | 273 | 69 | 202 | 205 | 283 | 69 | 208 | 205 | 289 |
| 27 | 179 | 68 | 199 | 204 | 281 | 71 | 207 | 212 | 291 | 71 | 213 | 212 | 297 |
| 28 | 186 | 70 | 203 | 211 | 289 | 73 | 211 | 219 | 299 | 73 | 217 | 219 | 305 |
| 29 | 193 | 71 | 207 | 218 | 297 | 74 | 215 | 226 | 307 | 74 | 221 | 226 | 313 |
| 30 | 199 | 72 | 212 | 225 | 305 | 75 | 220 | 233 | 315 | 75 | 226 | 233 | 321 |
| 31 | 206 | 73 | 217 | 234 | 315 | 76 | 225 | 242 | 325 | 76 | 231 | 242 | 331 |
| 32 | 213 | 75 | 221 | 243 | 323 | 78 | 229 | 251 | 333 | 78 | 235 | 251 | 339 |
| 33 | 219 | 77 | 225 | 252 | 331 | 80 | 233 | 260 | 341 | 80 | 239 | 260 | 347 |
| 34 | 226 | 78 | 229 | 261 | 339 | 81 | 237 | 269 | 349 | 81 | 243 | 269 | 355 |
| 35 | 232 | 80 | 233 | 268 | 347 | 83 | 241 | 278 | 357 | 83 | 247 | 278 | 363 |
| 36 | 239 | 82 | 237 | 279 | 355 | 85 | 245 | 287 | 365 | 85 | 251 | 287 | 371 |
| 37 | 246 | 83 | 240 | 288 | 365 | 86 | 248 | 296 | 375 | 86 | 254 | 296 | 381 |
| 38 | 252 | 85 | 244 | 297 | 375 | 88 | 252 | 305 | 385 | 88 | 258 | 305 | 391 |
| 39 | 259 | 86 | 248 | 306 | 385 | 89 | 256 | 314 | 395 | 89 | 262 | 314 | 401 |
| 40 | 266 | 88 | 252 | 315 | 395 | 91 | 260 | 323 | 405 | 91 | 266 | 323 | 411 |
| 41 | 272 | 90 | 255 | 324 | 405 | 93 | 263 | 332 | 415 | 93 | 269 | 332 | 421 |
| 42 | 279 | 91 | 258 | 333 | 415 | 94 | 266 | 341 | 425 | 94 | 272 | 341 | 431 |
| 43 | 286 | 93 | 261 | 342 | 425 | 96 | 269 | 350 | 435 | 96 | 275 | 350 | 441 |
| 44 | 292 | 95 | 264 | 351 | 435 | 98 | 272 | 359 | 445 | 98 | 278 | 359 | 451 |
| 45 | 299 | 96 | 268 | 360 | 444 | 99 | 275 | 368 | 454 | 99 | 281 | 368 | 460 |

Measurements for all sections are shown in mm.

To meet the clinical requirement of a patient Blatchford have the right to modify the UKCA marked footwear and by default this will be labeled as Custom-made device. Blatchford maintains a QMS system which is ISO 9001 and ISO 13485 accredited. Blatchford products and systems meet the essential requirements of EU Medical Devices Regulation 2017 (MDR).

Distributed by

Blatchford | KSA |  **Mendivil**

blatchfordmobility.com | orthoticscs@blatchford.co.uk | [@blatchfordgrp](https://twitter.com/blatchfordgrp)